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
APPAREIL A RAYONS X

Experiment topics:

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

SPECTRE ATOMIQUE (Atomic Spectra)


- Hydrogen Balmer Series
- Helium Spectrum
- Mercury Doublet

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The 3 brightest lines of the Balmer series for Hydrogen. A scan with the Light Sensor set on a higher sensitivity reveals a violet line which cannot be seen in this scan.

Advantage:

The open construction of this spectrophotometer allows the entire spectrum to be seen while the intensity versus angle is graphed in real-time. The wavelengths of the discrete lines of the atomic spectra of various gases are measured using a grating spectrophotometer. The atomic spectra of Hydrogen, Helium and Mercury are scanned by hand using a grating spectrophotometer, which measures relative light intensity as a function of angle. From the resulting graph, the wavelengths of the spectral lines are determined by measuring the angle from the central maximum to each line. First and second order lines are examined. The spectrum of Sodium (the Sodium doublet cannot be resolved) is used to calibrate the diffraction grating. The wavelengths of the spectral lines are compared to the accepted values and, in the case of Hydrogen, the electron orbit transitions corresponding to the lines are identified.

**LA CHARGE D'ELECTRON "MILLIKAN"
(Charge of an Electron).**

- Accurately Measure the Charge of a Single Electron
- Recreate Robert Millikan's Historical Experiment

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
Small droplets of oil are introduced into a chamber where an electric field of known strength is present. Using the viewing scope and a stopwatch, the velocity of a falling oil droplet is measured and recorded. Next, the electric field in the chamber is increased, causing the oil droplet to move upward. This allows the measurement of the force on the droplet and, ultimately, the charge of the droplet. By measuring the charge of several different oil droplets, the smallest difference in charge between them can be equated to the charge of an electron.

Advantage:


Charge of an Electron Experiment features a 30X, bright-field, erect-image microscope for clear viewing of the oil droplets. The droplet viewing chamber utilizes a special condenser to minimize droplet drift typically caused by droplet illumination and outside air currents. An ionization source allows the droplet charge to be changed.

Experiment Includes:

- Millikan Oil Drop Apparatus
- Basic Digital Multimeter
- High Voltage Power Supply
- Large Rod Base
- 45 cm Steel Rod (2)
- Banana Plug Cord-Red (5 Pack)
- Banana Plug Cord-Black (5 Pack)
- Charge of an Electron Experiment Manual

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Bench for experiments with radioactivity:

Bench for Radioactive Sources and GM-tubes, providing a stable and convenient setup, allowing reproducible results when making experiments with radioactivity. Investigate how the radioactive radiation spreads from the source in both angle and distance or investigate how radiation is absorbed by lead and aluminum.

Comprises a 60 cm long bench, with a mm-scale, holder for GM-tubes and holder for radioactive Sources. The later can be fitted with a diameter 6 mm steel rod, for hanging absorber plates in front of the radioactive source. The included absorber plates are supplied in a convenient storage box, and includes: 5 ea. 2mm and 10 ea. 1 mm thick lead plates plus 2 ea. 3 mm; 6 ea. 2mm and 4 ea. 0,5 mm Aluminum plates. The lead-plates have been enameled to protect fingers from lead contamination. The holder for the radioactive sources can be turned, and is fitted with a degree scale, +/- 110 degree, 1 degree resolution.

DEVIATION DE LA PARTICULE BETA

Deflection of beta particles:

The deflection of beta particles is easily demonstrated with this accessory for the bench for experiments with radioactivity (5141.00). The magnetic field is provided by a pair of strong permanent magnets. The direction of the field can be reversed and the magnets can be removed altogether.

A beta source is mounted in a collimator that can turn around the centre of the magnets, thereby changing the angle between the incoming particle beam and a Geiger-Müller detector.

To demonstrate the effect of the magnetic field on the beta radiation, you start with the magnets removed from the apparatus and align the beta source and the detector. When the magnets are re-introduced, the radiation “disappears”. By changing the angle between source and deflector, the radiation can be found again – but only if the magnetic field is oriented correctly.

From the angle of deflection the kinetic energy of the particles can be deduced. This way the energy spectrum can be sketched to show its continuous nature.